


Physical Limitations for Electron Microscopy of Lithium and Sulfur Battery Materials

Barnaby D.A. Levin, David A. Muller
Cornell University Applied Physics

Why are Lithium and Sulfur difficult to image in TEM? How can these difficulties be overcome?


The Electron Beam Will Damage Lithium

- There is a cross section for radiation damage.
- There is also a cross section for imaging signal.
- These cross sections depend on electron beam energy
- Our goal is to maximize the number of imaging events per damage event.



- References**
1. J. R. Breedlove, J.; Trammell, G. T., *Science* 1970, 170, 1310-1313.
 2. Henderson, R. *Quarterly Reviews of Biophysics* 1995, 28, 171-193.
 3. Bradley, C. R. *Calculations of Atomic Sputtering and Displacement Cross-Sections in Solid Elements by Electrons with Energies from Threshold to 1.5MV*; Argonne National Lab: 1988.
 4. Egerton, R. F., *Electron Energy Loss Spectroscopy in the Electron Microscope*. 3rd ed.; Springer: New York, 2011.

Two principal forms of electron beam induced damage.



Damage Cross Section (σ_D) for Li peaks at low beam energy.³
Lower bound (Mckinley-Fesbach) assumes no inter-atom interaction. Upper bound (Kinchin-Pease) assumes cascades induced.

Signal cross sections are larger at lower beam energy.⁴
Cross sections for K-edge (σ_K), Plasmon (σ_P) and elastic scattering (σ_E) plotted with σ_D .


Ratios imply imaging conditions for Li are superior at lower beam energy.
Ratios of σ_K and σ_E against σ_D are greater at lower voltages, yielding more imaging events per damage event.

Conclusions

- Imaging with Low Voltage TEM can reduce knock-on damage to Li.
- Voltage also subject to other factors e.g. sample thickness, beam spreading.

Bulk Sulfur Will Sublimate in the Vacuum of a TEM

Sulfur vapor pressure is $\sim 10^{-6}$ Torr at room temperature.⁵
Typical TEM pressures are $\sim 10^{-8}$ Torr.




Vapor pressure curve for Sulfur, Kelvin/Antoine formula⁶ extrapolation from experimental data. Typical TEM sample pressures $\sim 10^{-8}$ Torr, (red line).

There are 3 ways to avoid Sulfur sublimation:

1. Encapsulate in material with $<$ a few nm sized pores.¹⁰
2. Use environmental or air electron microscope.
3. Use a Cryo-electron microscope.


Solution 1: Encapsulation



HAADF STEM image of Sulfur encapsulated in a PAN shell.
Inset: EDX map. Red = Carbon, Green = Sulfur.
Sample provided by: Yingchao Yu, Weidong Zhou, Cornell Chemistry

Solution 2: Image in Air

Going from low-mag 200keV TEM to 30keV Air STEM, there is no significant loss of resolution.




a ABF STEM image of Sulfur particles in PAN shell.
b The sample was transferred to Air-SEM and imaged again with a STEM detector. Resolution sufficient to resolve empty shells clearly ~ 10 nm.

Air SEM Imaging assistance from Kayla Nguyen, Muller Group, Cornell.

Conclusions



- At low magnification, Air-STEM can be used to image Sulfur particles without sublimation. Further improvements possible.
- For high magnification imaging, Cryo-electron microscopy is an option.

Acknowledgements for assistance and useful discussions, Muller Group, Cornell University.



References

5. Ferreira, a. G. M., & Lobo, L. Q. (2011). *J. Chem. Therm.* 43(2), 95–104.
6. Blundell & Blundell (2006), *Concepts in Thermal Physics*, OUP.
7. C. Briske, N.H. Hartshorne, D.R. Stranks, (1960) *J. Chem. Soc.* 1200–1209.
8. K. Neumann, (1934) *Z. Phys. Chem.* 171 416–420
9. G. Foutetier, (1944) *Compt. Rend.* 218 194–196
10. E.P. Barrett, L.G. Joyner, P.H. Halenda, (1951) *J. Am. Chem. Soc.*, 73 373-380



We would like to thank John Grazul and Mick Thomas for assistance in the TEM facility. Supported by the Energy Materials Center at Cornell, an Energy Frontier Research Center funded by the U.S. Department of Energy, Office of Science, BES Award DE-SC0001086. EM Facility support from the NSF MRSEC program (DMR 1120296).



Cornell University
College of Engineering
School of Applied and Engineering Physics